

## Quick Guide

1 Download the IBEX app.



2 Launch the IBEX app and turn on the probe.

3 Connect the probe with the smart device.



Get help from us



## E.I. MEDICAL IMAGING

E.I. Medical Imaging | [www.eimedical.com](http://www.eimedical.com)

[Info@eimedical.com](mailto:Info@eimedical.com)

815 14th St. SW, Unit C210

Loveland, CO 80537

# IBEX Ultrasound Imaging System

## Quick Guide

E.I. Medical Imaging | [www.eimedical.com](http://www.eimedical.com)  
Info@eimedical.com  
815 14th St. SW, Unit C210  
Loveland, CO 80537

1 Download the IBEX Ultrasound app.



Launch the IBEX Ultrasound app and turn on the probe.



3 Long press the power button until the light turns purple (about 3 seconds) to power on.



Connect the probe to the device via Wi-Fi. For the first time, scan the QR code on the probe to connect it.



Get help from us



# Scanning Basics

