

THE NEW LITENXT



IT'S HERE...THE NEW IBEX LITENXT

For 40 years, we've manufactured animal-specific ultrasound right here in the US. Now we've done it again..the next generation of bovine ultrasound is here!

The **LITENXT** sets a new standard with image quality never seen before in a bovine ultrasound.

- 30% smaller than traditional LITE ultrasound
- New **HYBRID eCli6** broadband probe —hardwired, combines linear and curved linear into one probe
- New **InSite NXT** headset with latest in display technology—OLEDs with 2x the resolution & 3x the brightness...70% larger image size
- Faster frame rate—up to 60fps
- Built-in wireless streaming to your phone or tablet
- Auto stand-by mode—allows 20% savings when not scanning, 5 hours continuous scanning on 1 battery
- Completely redesigned image software with the latest technology—tissue harmonics, synthetic imaging, speckle reduction—simplified bovine-specific exam presets
- Legendary IBEX ruggedness—designed and tested for the harshest field environment; designed, manufactured and supported in Colorado

E.I. MEDICAL IMAGING
1.866.365.6596 • www.eimedical.com

eCLi6 HYBRID TRANSDUCER



The new, patent-pending **eCLi6** transducer is hardwired to our latest **IBEX LITENXT** ultrasound machine and is a hybrid probe. The **eCLi6** transducer is unlike anything you've seen before!

This unique rectal probe offers the resolution of a flat linear with the deeper penetration and wider field-of-view of a convex model. The short array and tight pitch optimize image quality. Fits in your hand like a linear probe, while the very slight curvature makes ideal contact with any surface in a minimal-height form factor.

Based on the overwhelmingly positive response since we started testing it last fall, this transducer is a game-changer for both bovine and equine reproduction. Outperforms any standard rectal probe for image clarity, arm-in or arm-free.

See for yourself—free demo, no obligation! Call us at 1.866.365.6596 or go to www.eimedical.com/litenxt...you'll be glad you did!



TESTIMONIAL

"The LITENXT is so much clearer and better resolution than my IBEX PRO was. The small size makes it easy to just stick it in your pocket and go as opposed to having to wear a unit around your waist. The WiFi connection to your phone is a great option for letting the client see what I am seeing and even snapping a screenshot of the image on your phone to text to someone." *DR. PHIL LOWE*

LITENXT Rev.2.1

A WORD ABOUT OUR CABLES...

IBEX® cables tested to 1.2 million flex cycles without breaking

- Custom designed cables are abrasion resistant
- Only manufacturer to design all transducers specifically for animal use
- Kevlar reinforced



Made in Colorado, USA



STREAMING NOW BUILT-IN...

E.I. MEDICAL IMAGING

E.I. Medical Imaging® was founded in 1984 with the vision of developing state-of-the-art, real-time ultrasound to serve veterinarians and livestock producers world-wide.

Our products evolved with the needs of the market. We are the only manufacturer of portable ultrasound systems engineered and serviced in the United States specifically for the animal industry.

E.I. Medical Imaging—ultrasound solutions for your practice. Made in the USA since 1984.

CONTACT US

1.866.365.6596 [toll free]
1.970.669.1793
1.970.669.1902 [fax]
info@eimedical.com
www.eimedical.com

