

THE NEW PRONXT



IT'S HERE...THE ALL-NEW PRONXT

For 40 years, we've manufactured animal-specific ultrasound right here in the US. Now we've done it again..the next generation of portable ultrasound is here!

The **PRONXT** sets a new standard with image quality never seen before in a large animal—bovine and equine—repro ultrasound.

- 25% larger, high-resolution sunlight readable display
- Simplified, sealed keyboard for ease of use and cleaning
- New **HYBRID eCli6** broadband probe —hardwired, combines linear and curved linear into one probe
- Faster frame rate—up to 60fps
- Built-in wireless streaming to your phone or tablet
- Auto stand-by mode—allows 20% savings when not scanning, 4 hours continuous scanning on 1 battery
- Completely redesigned image software with the latest technology—tissue harmonics, synthetic imaging, speckle reduction—significant image quality improvement
- Compatible with InSite NXT goggles [optional]
- Legendary IBEX ruggedness—designed and tested for the harshest field environment; designed, manufactured and supported in Colorado

E.I. MEDICAL IMAGING
1.866.365.6596 • www.eimedical.com

eCLi6 HYBRID TRANSDUCER



The new, patent-pending **eCLi6** transducer is hardwired to our latest **IBEX PRONXT** ultrasound machine and is a hybrid probe. The **eCLi6** transducer is unlike anything you've seen before!

This unique rectal probe offers the resolution of a flat linear with the deeper penetration and wider field-of-view of a convex model. The short array and tight pitch optimize image quality, while the very slight curvature makes ideal contact with the mucosal surface in a minimal-height form factor.

Based on the overwhelmingly positive response since we started testing it last fall, this transducer is a game-changer for bovine reproduction. Outperforms any standard rectal probe for image clarity, arm-in or arm-free.

See for yourself—free demo, no obligation! Call us at 1.866.365.6596 or go to www.eimedical.com/litenxt...you'll be glad you did!



TESTIMONIAL

Fantastic upgrade from the PRO series aimed at improving the quality of the services I provide to cow/calf clients. Image quality is cutting edge, keyboard is very simple and intuitive, all at a price that works for cattle veterinarians. *TK Livestock Services*

A WORD ABOUT OUR CABLES...

IBEX® cables tested to 1.2 million flex cycles without breaking

- Custom designed cables are abrasion resistant
- Only manufacturer to design all transducers specifically for animal use
- Kevlar reinforced



E.I. MEDICAL IMAGING

E.I. Medical Imaging® was founded in 1984 with the vision of developing state-of-the-art, real-time ultrasound to serve veterinarians and livestock producers world-wide.

Our products evolved with the needs of the market. We are the only manufacturer of portable ultrasound systems engineered and serviced in the United States specifically for the animal industry.

E.I. Medical Imaging—ultrasound solutions for your practice. Made in the USA since 1984.

CONTACT US

1.866.365.6596 [toll free]

1.970.669.1793

1.970.669.1902 [fax]

info@eimedical.com

www.eimedical.com

