

THE NEW SA²



IT'S HERE...THE NEW SA²

For 40 years, we've manufactured animal-specific ultrasound right here in the US. Now we've done it again... the next generation of the SA².

With new broadband probes, improved imaging and battery performance, the SA² provides high-quality ultrasound in a rugged, portable machine. Great for the clinic or field!

- New broadband probes for all species
- Completely redesigned image software with the latest technology
- Customized, veterinary-specific exam presets
- Sunlight-readable LED display—44% larger XGA screen
- WiFi streaming via EVOStream app
- Ruggedized, splash resistant, and portable [just over 7 lbs]—sealed keyboard and probe connections allow easy cleaning
- Interchangeable batteries for continuous use
- Designed and manufactured in Colorado

E.I. MEDICAL IMAGING
1.866.365.6596 • www.eimedical.com/sa2



BROADBAND TRANSDUCERS



eCL3
2-5 MHz macro convex



eCLi4
3-7 MHz curved rectal



eCLi6 NEW!
3-9.5 MHz hybrid curved linear



eC6
4-9 MHz micro convex



eC9
5-10 MHz micro convex



eL7
5-9 MHz linear rrectal



eL14
6-14 MHz tendon



eC90PU
5-10 MHz broadband array for transvaginal oocyte aspiration in horses and cattle. Adjustable handle, small height, lightweight.



A WORD ABOUT OUR CABLES...

IBEX® cables tested to 1.2 million flex cycles without breaking

- Custom designed cables are abrasion resistant
- Only manufacturer to design all transducers specifically for animal use
- Kevlar reinforced

...AND OUR NEW TRANSDUCERS

As part of our commitment to innovation, we have re-designed our probes for improved penetration and resolution. Broadband, variable frequency, our transducers are designed specifically for rugged environments. But if the worst happens, they can be repaired right here in our Colorado facility.

E.I. MEDICAL IMAGING

E.I. Medical Imaging® was founded in 1984 with the vision of developing state-of-the-art, real-time ultrasound to serve veterinarians and livestock producers world-wide.

Our products evolved with the needs of the market. We are the only manufacturer of portable ultrasound systems engineered and serviced in the United States specifically for the animal industry.

CONTACT US

1.866.365.6596 [toll free]
1.970.669.1793
1.970.669.1902 [fax]
info@eimedical.com
www.eimedical.com



Made in Colorado, USA

