


i3 GOGGLES



The all-new *i3* goggles* are offered in 2 models... monocular and binocular.

Both versions come with a pre-attached hat clip, designed to be clipped to the bill of a hat or glasses.

A head strap is also included with the binocular version, along with hardware for mounting.

All models are include:


- 2 sets of rubber glare shields, one with a flare to shield sunlight and one without for better peripheral vision
- User-replaceable cables

Both models are small and lightweight, and are compatible with Ibex EVO, Pro, and Lite as well as with Bantam and Easi Scan 1-4.



Custom designed optics

- 800x600 color Oled display(s)
- Half the weight of InSight 2 goggles
- Specially designed eye cups for inside and outdoor viewing
- Adjustable links, universal clip
- GoPro mount
- Standard, extended warranties available [regular use]



For more information or to schedule a demonstration,
contact us at 866.365.6596 • info@eimedical.com

**patent pending*

E.I. Medical Imaging®